


Size Exclusion Chromatography (SEC)


PolyHYDROXYETHYL A™ can be used for SEC from the largest proteins to amino acids. It is noteworthy for the dramatic shift in fractionation range with a mildly denaturing mobile phase such as 50 mM formic acid.



北京金欧亚科技发展有限公司

北京崇文区左安门内大街8号伟图大厦301室 PC: 100061
Tel: 010-67136152/67100708 Fax: 010-67114016/67113925
<http://www.chromatogr.com> E-mail: china.hplc@163.com

The chaotrope in the mobile phase disrupts the hydrogen bonding between the polymer chains in the coating, increasing their steric radius and making the space between them permeable. This permits the fractionation of solutes as small as dipeptides and amino acids.




Our columns have been used to assess the size of the small peptides present in commercial **protein hydrolyzates** as well as the products from **proteasome** digestion.

Our columns work well with proteins.

北京金欧亚科技发展有限公司

北京崇文区左安门内大街 8 号伟图大厦 301 室 PC: 100061
Tel: 010-67136152/67100708 Fax: 010-67114016/67113925
<http://www.chromatogr.com> E-mail: china.hplc@163.com

PolyHYDROXYETHYL A, 1000-Å
SEC with 50 mM Formic Acid



Flow rate and column size: Optimum resolution is obtained with a 200x9.4-mm column at a flow rate around 0.5-0.6 ml/min. Since most HPLC systems can deliver such flow rates reliably, this is our most popular size SEC column. Optimum resolution with a 200x4.6-mm column would be obtained with a flow rate of 130-140 μ l/min. Make sure that your HPLC system can deliver such flow rates reliably, and minimize dead volume in the connections and detector flow cell (use a microbore-scale cell).

北京金欧亚科技发展有限公司

北京崇文区左安门内大街8号伟图大厦301室 PC: 100061
 Tel: 010-67136152/67100708 Fax: 010-67114016/67113925
<http://www.chromatogr.com> E-mail: china.hplc@163.com